



KREO HMI TUTORIAL
FDA – CFR21

Tutorial dedicated to the development of an application FDA
cfr21 part 11 compliant

Connect
Ideas.
Shape
solutions.



Introduction

This document provides a brief description of FDA CFR21 part 11 and the procedures in order to develop a KREO HMI application compliant with this regulation.

The purpose is to have the legal equivalence of electronic documents (digital records and electronic signature) compared to traditional paper ones.

In order to be compliant with the CFR21 Part 11 standard, it is necessary to ensure that the recorded data always refer to the operator (Electronic Signature) and, in addition, specific policies are implemented that make it impossible to modify the electronic recorded data.

The regulation

FDA (Food and Drug Administration) is the U.S. agency responsible for controlling and regulating production processes in the food, pharmaceutical, and chemical industries.

It is not only American companies that are subject to its control, but all those that operate in the US and export there.

CFR 21 is a regulation issued by the FDA in 1997 in collaboration with the U.S. government regarding the use of technology in specific manufacturing process procedures. Part 11 is divided into two main sections:

- Electronic Archives
- Electronic Signatures

CFR 21 Part 11 regulates the way in which data stored in electronic support is handled and the related security issues.


The objective of CFR 21 Part 11 is therefore to verify that all procedures adopted in the process comply with the requirements.

CFR 21 Part 11 of the FDA is also a guarantee tool for companies such as ESA SPA, which want to offer their customers high and consistent quality and safety standards in their products.




How to do:

- 1) In order to enable the FDA functionalities it is necessary to open the Configuration – FDA section, select the Enable option and then insert the general settings:





3) Now considering a typical RUNTIME page, the interaction with the different objects to be FDA recorded (numeric field, login, logout, recipe download, language change,...) will be displayed the log below:



FDA

Index	Time	Code	Object	User	FDA Note	TagCurrentValue	TagNewValue	ElectronicSign	FDA Comment
	17:24:29 30/6/2021	System / Startup							
1	17:24:34 30/6/2021	Users / Logout		defaultuser					
2	17:24:34 30/6/2021	Users / Login		defaultuser					
3	17:25:13 30/6/2021	UI / Input confirmed	w470607	defaultuser	ChangingTag	0	23	Adam Smith	change volume
4	17:25:33 30/6/2021	UI / Button pressed	w5306011	defaultuser	LoginAccess			Admin1	log1
5	17:25:33 30/6/2021	Users / Login		Admin					
6	17:26:35 30/6/2021	UI / Button pressed	w53034017	Admin	Recipe download			P. Rossi	rec1
7	17:26:59 30/6/2021	UI / Input entered	w53050015	Admin	ChangeLanguage				lang
8	17:27:18 30/6/2021	UI / Button pressed	w5207504	Admin	ExportingFDA			me	exp1

23 Volume (litres)





4) At each click on the field/key/image object you will get the input request for description and signature (if enabled in the project events):

The screenshot shows a web browser window displaying an FDA interface. The main content is a table with the following columns: Index, Time, Code, Object, User, FDA Note, TagCurrentValue, TagNewValue, ElectronicSign, and FDA Comment. The table contains several rows of data, with the last row (Index 8) highlighted in blue. A modal dialog box is overlaid on the table, containing two input fields: "Description" and "Electronic Signature". Below the input fields are two buttons: "Cancel" and "OK".

Index	Time	Code	Object	User	FDA Note	TagCurrentValue	TagNewValue	ElectronicSign	FDA Comment
	17:24:29 30/6/2021	System / Startup							
1	17:24:34 30/6/2021	Users / Logout		defaultuser					
2	17:24:34 30/6/2021	Users / Login		defaultuser					
3	17:25:13 30/6/2021	UI / Input confirmed	w470607	defaultuser	ChangingTag	0	23	Adam Smith	change volume
4	17:25:33 30/6/2021	UI / Button pressed	w					Admin1	log1
5	17:25:33 30/6/2021	Users / Login							
6	17:26:35 30/6/2021	UI / Button pressed	w					P. Rossi	rec1
7	17:26:59 30/6/2021	UI / Input entered	w						lang
8	17:27:18 30/6/2021	UI / Button pressed	w					me	exp1

23 Volume (litres)

login logout Send Recipe  export FDA



5) The same log can be exported in PDF / CSV format through the function: FDATracingExport. Below is an example of the result:

FDA.pdf - Adobe Acrobat Reader DC (32-bit)
File Modifica Vista Firma Finestra Aiuto

Home Strumenti FDA.pdf x

File Type **RUNTIME_EVENTS_EXPORT**

Station Name **EW410N_MAXWELL**

Project Id **testMANUALE**

FW Version **200**

Creation Date **30-06-2021**

Creation Time **17:27:18**

Id	Time User	Module/Action Signature	Object Note	Comment
1	30-06-2021 17:24:29.854	System - Startup		
2	30-06-2021 17:24:34.921 defaultuser	Users - Logout		
3	30-06-2021 17:24:34.922 defaultuser	Users - Login		
4	30-06-2021 17:25:13.754 defaultuser	UI - Input confirmed Adam Smith	w470607 ChangingTag	change volume
5	30-06-2021 17:25:33.086 defaultuser	UI - Button pressed Admin1	w5306011 LoginAccess	log1
6	30-06-2021 17:25:33.159 Admin	Users - Login		
7	30-06-2021 17:26:35.973 Admin	UI - Button pressed P. Rossi	w53034017 Recipe download	rec1
8	30-06-2021 17:26:59.260 Admin	UI - Input entered	w53050015 ChangeLanguage	lang
9	30-06-2021 17:27:18.293 Admin	UI - Button pressed me	w5207504 ExportingFDA	exp1



The FDA events:

Events		>>
OnStart	None	
OnFdaError	None	
OnActivityOn	None	
OnActivityOff	None	
OnStop	None	
OnError	None	
OnAnyUserLogin	None	
OnAnyUserLogout	None	
OnAnyUserLoginError	None	
OnAnyUserInfoChanged	None	
OnAnyUserCreated	None	
OnAnyUserDeleted	None	
OnAnyUserLocked	None	
OnAnyUserUnlocked	None	
OnUsersReset	None	



Connect
ideas.
shape
solutions.

[ESA S.p.A. | www.esa-automation.com](http://www.esa-automation.com) |